

Sony Cell Sorter MA900

Default filter setting

Detector	Filter	Fluorochrome
FL1	525/50	FITC, BB515, AF488, GFP, YFP, mCitrine, CFSE
FL2	585/30	PE,dsRed,ddTomato
FL3	617/30	mCherry, PE-TR, PE-Dazzle594, PI
FL4	695/40	mPlum,7-AAD,PeCy5,PerCP,Pe-Cy5.5,PerCP-Cy5.5, PerCP-eF710
FL5	785/60	Pe-Cy7
FL6	450/50	BV421, Pacific Blue, DAPI, AF405, mCFP,Hoechst33342
FL7	525/50	AmCyan, BV510
FL8	585/30	BV570
FL9	617/30	BV605
FL10	665/30	BV650, APC, AF647
FL11	720/60	BV711,APC-Cy5.5,AF700
FL12	785/60	BV785,APC-Cy7,APC-AF750

